

Baaritaanka Rapid Antigen Test (RAT)

Sida loogu isticmaalo baaritaanka rapid antigen test si sax ah.

Dhamaan baaritaanada rapid antigen tests, waxaa waajiba kugu ah inaad mar walba:

- Hubiso baaritaanku in aanu dhicin
- Dhaq ama nadiifi gacmahaaga ka hor iyo ka dib markaad qaadayso baaritaanka
- Si degan u raac awaamiirta baaritaanka ee lagu muujiyey awaamiirta ka muuqata baakadaynta



Waxay tahay inaad samayso ama aadan samaynayn



Baaritaanada isticmaala candhuufada ama calyada

30 daqiiqo ka hor inta aadan qaadin baaritaanka, ha:



Cunin,



Dhuuqin/
cabin



Ilkahaaga
ha cadaysan



Calaashan
Xanjo



Ka hor inta aadan qaadin baaritaanka, ku qufac mid xoogan 3 ilaa 5 jeer maaskiga ama tii-shiyuu kana fogow dadka.



Baaritaanada isticmaala sambalada sankaa



Duufso ka hor inta aadan samayn baaritaanka.



Ha taaban dhinaca jilicsan oo qoriga yar, kaas oo loo yaqaan swab.

Xasuuso, baaritaan walba waxaa la isticmaali karaa keliya hal mar.



Macluumaad intaas ka badan oo ku saabsan Rapid Antigen Testing, wac

National Coronavirus Helpline: 1800 020 080

Dooro xalka 8aad caawimada turjumaan lacag la'aana.

Waxaad kaloo booqan kartaa health.gov.au



Australian Government